








Product name:	OD Eleva All Digital (OmniDiagnost Eleva All Digital)	
Identification code	708-028	
Total weight (in Kg)	2495 kg (approximately; dependent on specific configuration)	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	Veenpluis 6
	Zip code:	5684 PC Best
	Country:	Netherlands
	Electronic info:	http://www.healthcare.philips.com/us/about/sustainability/recycling/

Recycle Info	Items:	Location
Special attention 	<ul style="list-style-type: none"> Be aware of possibly contaminated system parts and materials! (biological hazard) For dismantling activities Treatment Facilities must consider the national requirements. For personnel that can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account 	System parts that were in the patient environment, and that were not disinfected
	<ul style="list-style-type: none"> Removal of units / weights can cause the system(parts) to tilt! 	
	<ul style="list-style-type: none"> Removal of units / weights can cause unexpected movements of guidances! 	
	<ul style="list-style-type: none"> Release of brakes can cause unexpected movements of guidances! Brakes cannot prevent unexpected movements due to the removal of units /weights! 	
	<ul style="list-style-type: none"> High-voltage parts (e.g. capacitors) are marked with  	
	<ul style="list-style-type: none"> Before dismantling the vacuum II-Insert, drill a small hole to let air in the insert 	II-TV (page 4-12)
	<ul style="list-style-type: none"> Vacuum glass tube can implode! 	X-ray tube (page 14-16)
	<ul style="list-style-type: none"> When present: take caution dismantling a CRT screen 	CRT screens (page 27-36)
Fluids / Gases 	Items: <ul style="list-style-type: none"> Cooling liquidGlycoshell (in cooling unit, when present) [not for air-cooled X-ray tube versions] 	
	<ul style="list-style-type: none"> Transformer oil, type: Shell Diala 	Cabinets (page 37-43)
	<ul style="list-style-type: none"> Oil 	X-ray tube (page 14-16)


Batteries  To be Removed	Type:	Location
	Battery, 4x alkaline 1,5V [44 grams] (per “remote control/viewpad”)	 (page 3)
	CR2032 3.0V Lithium coin cell of 3.2 gram (“Dell PC”)	Electronics (page 16)
	Alkaline battery	Cabinets (page 37-43)
	1x CR2032 3.0V Lithium coin cell	Cabinets (page 37-43)
	1x 3.5V Lithium battery	Cabinets (page 37-43)
	Lithium chrome cell, 3V (2x)	Cabinets (page 37-43)
Hazardous  To be Removed	Substances:	Location
	Lead (Pb) for X-ray shielding	II-TV (page 4-12) Grid (page 13) X-ray tube (page 14-16) Collimator (page 17)
	Lead (Pb) for soldering	Electronics (page 18) Display screens (page 19-36) Cabinets (page 37-43)
	Cadmium (Cd) + Beryllium Oxide (BeO) inside the II-Insert on the glass output window	II-TV (page 4-12)
	Beryllium Copper (BeCu)	Electronics (page 18)
	Beryllium Copper (BeCu) contact springs	Cabinets (page 37-43)
	Mercury (Hg) in switch on printed circuit board for systems delivered before September 2006	II-TV (page 4-12)
	Mercury (Hg) in specific LCD screens, when these LCD screens are present	LCD screens (page 19-26)

Note: to facilitate recycling, all plastic parts weighing > 50 grams are marked according to ISO11469 & ISO1043.

Display screens

OmniDiagnost Eleva All Digital


When present:



X-ray tube

Collimator

One or more remote controls /viewpads:



II-TV

Grid (inside)

Electronics (inside)

Cabinets

II-TV: page 4-12

Grid: page 13

X-ray tube: page 14-16

Collimator: page 17

Electronics: page 18



Display screens: page 19-36


Cabinets: page 37-43

Title: Recycling passport OD Eleva All Digital 708-028

DocID: XDR054-090790


23cm II-TV:


Recycle Info	Items:	Location
Special attention 	Before dismantling the vacuum II-Insert, drill a small hole to let air in the insert	
Hazardous	Substances:	Location
 To be Removed	Lead (Pb)	1, page 4-5 + 7
	Cadmium (Cd) + Beryllium Oxide (BeO) inside the II-Insert on the glass output window	2, page 4-5 + 7
	Beryllium Copper (BeCu) contact springs	3, page 6
	Mercury (Hg) in switch on printed circuit board for systems delivered before September 2006	4, page 4 -5 + 8
	Pb is present in the soldering process of some PCBs	




Before September 2006:

See detail B
Printed circuit board







IDSC RACK 4522 163 2455X in cabinet



Detail A






Detail B printed circuit board 4522 167 02681 up and including 4522 167 02687




Detail B printed circuit board 4522 167 02431 up and including 4522 167 02439


31cm & 38cm II-TV:

Recycle Info	Items:	Location
Special attention 	Before dismantling the vacuum II-Insert, drill a small hole to let air in the insert	
Hazardous	Substances:	Location
 To be Removed	Lead (Pb)	1, page 9-10
	Cadmium (Cd) + Beryllium Oxide (BeO) inside the II-Insert on the glass output window	2, page 9-10
	Beryllium Copper (BeCu) contact springs	3, page 11
	Mercury (Hg) in switch on printed circuit board for systems delivered before September 2006	4, page 9 + 12
	Pb is present in the soldering process of some PCBs	




Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790






IDSC RACK 4522 163 2455X in cabinet





Detail B printed circuit board 4522 167 03471 up and including 4522 167 03475



Detail B printed circuit board 4522 167 02431 up and including 4522 167 02439

Grid:

Hazardous  To be Removed	Substances:	Location
	Lead (Pb 99,5%)	Enclosed between cover plates






Cross-section of grid:



**Paper
Interspace
Material**


Lead

X-ray tube SRM 0608 ROT-GS505 / 9890-000-8518x

Recycle Info	Items:	Location
Special attention 	Vacuum glass tube can implode!	X-ray tube
Fluids / Gases	Items:	Location
	Oil: 4,2 kg HV-transformer oil contains no PCBs	
Hazardous	Substances:	Location
	Lead as x-ray shielding inside housing: 5,5 kg	Figure below
To be Removed	Beryllium as x-ray window in vacuum envelope: 1 g	Figure below

Beryllium as x-ray window
in vacuum envelope




Lead as x-ray shielding
inside housing




Vacuum glass tube can implode!

steel, iron	iron, low alloy (<5%)	3,1KG
	iron, high alloy (>5%)	0,7KG
steel, iron		3,8KG
nonferrous metals and alloys	aluminium and -alloys	6,8KG
	copper and -alloys	2,7KG
	molybdenum and -alloys	0,9KG
nonferrous metals and alloys		10,4KG
glass / ceramics	ferrite	0,3KG
	glass	0,1KG
glass / ceramics		0,4KG
plastics / organic substances	oil	4,2KG
	thermoset	2KG
	elastomer	0,1KG
plastics / organic substances		6,3KG
relevant materials	lead and -alloys	5,5KG
	tungsten and -alloys	0,1KG
	beryllium and -alloys	0,001KG
relevant materials		5,601KG

X-ray tube SRO33100-R.351,0-3/90,150 / 9890-000-8585x




Recycle Info	Items:	Location
Special attention 	Vacuum glass tube can implode!	X-ray tube
Fluids / Gases	Items:	Location
	Oil: 3,5 kg HV-transformer oil contains no PCBs	
Hazardous	Substances:	Location
 To be Removed	Lead as x-ray shielding inside housing: 5,1 kg	Figure below




steel, iron	iron, low alloy (<5%)	2,86KG
	iron, high alloy (>5%)	0,2KG
steel, iron		3,06KG
nonferrous metals and alloys	aluminium and -alloys	4,3KG
	copper and -alloys	2,2KG
	molybdenum and -alloys	0,5KG
nonferrous metals and alloys		7KG
glass / ceramics	glass	0,2KG
glass / ceramics		0,2KG
plastics / organic substances	oil	3,5KG
	thermoset	1,9KG
	elastomer	0,1KG
plastics / organic substances		5,5KG
relevant materials	lead and -alloys	5,1KG
	cobalt and -alloys	0,24KG
	tungsten and -alloys	0,08KG
relevant materials		5,42KG

Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790

X-ray tube SRO33100-R.360,0-3/90,150 / 9890-000-8610x

Recycle Info	Items:	Location
Special attention 	Vacuum glass tube can implode!	X-ray tube
Fluids / Gases	Items:	Location
	Oil: 3,8 kg HV-transformer oil contains no PCBs	
Hazardous	Substances:	Location
 To be Removed	Lead as x-ray shielding inside housing: 3,8 kg	Figure below






Vacuum glass tube can implode!

steel, iron	iron, high alloy (>5%)	4,3KG
	iron, low alloy (<5%)	3,1KG
steel, iron		7,4KG
nonferrous metals and alloys	copper and -alloys	2,2KG
	aluminium and -alloys	1,7KG
	molybdenum and -alloys	0,5KG
nonferrous metals and alloys		4,4KG
glass / ceramics	glass	0,2KG
glass / ceramics		0,2KG
plastics / organic substances	oil	3,8KG
	thermoset	1,4KG
	elastomer	0,1KG
plastics / organic substances		5,3KG
relevant materials	lead and and -alloys	3,8KG
	tungsten and -alloys	0,08KG
relevant materials		3,88KG



Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790

Collimator:

Hazardous  To be Removed	Substances: Lead (Pb 99,5%); 3.7 – 4.1 kg; glued to several parts in the collimator	Location 1; see next pictures
---	---	---




Electronics:

Batteries  To be Removed	Type:	Location
	CR2032 3.0V Lithium coin cell of 3.2 gram ("Dell PersonalComputer" [Philips-indication: Viewforum hardware])	Inside Dell PC
Hazardous  To be Removed	Substances:	Location
	BerylliumCopper (BeCu)	Contact springs between hardware-racks; see also cabinets (page 38-44)
	Lead (Pb) is present in the soldering of some PCBs	PCBs (PrintedCircuitBoards)



Display screens:

FOLLOWING PAGES PROVIDE INFORMATION ON VARIOUS SCREENS POSSIBLY PRESENT IN THE SYSTEM.

LCD screen FIMI MML1801-IP1P / 9919-320-5067x | PAGE 1 of 2

Hazardous  To be Removed	Substances:		Location
	Type	Quantity	
	Cd	0	
	Hg	21 mg max. (Mercury is present in the backlight lamps)	Next figure (1)
	Pb	Lead is present in the soldering process of PCBs	Next figure (2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	


LCD screen FIMI MML1801-IP1P / 9919-320-5067x | PAGE 2 of 2




Material		
Fe	6 kg	-
Al	0	-
Cu	0.11 kg	Cables
Plastics	1.5 kg	-
Boards (S ² > 10cm ²)	96 cm ² / 260 g 320 cm ² / 230 g 144 cm ² / 67 g 100 cm ² / 50 g 46 cm ² / 40 g	S.M.P.S. (item 2.1 in the picture) Logic Board (item 2.2 in the picture) Inverter (item 2.3 in the picture) LCD Driver (item 2.4 in the picture) PCB Light Sensor (item 2.5 in the picture)
LCD	3.97 kg	18"

Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790

LCD screen FIMI MML1802-IP10 / 9919-320-5069x


Hazardous  To be Removed	Substances:		Location
	Type	Quantity	
	Cd	0	
	Hg	21 mg max. (Mercury is present in the backlight lamps)	Next figure (1)
	Pb	Lead is present in the soldering process of PCBs	Next figure (2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	





Material		
Fe	2.3 kg	-
Al	0	-
Cu	0.1 kg	Cables
Plastics	1 kg	-
Boards (S ² > 10cm ²)	96 cm ² / 260 g 320 cm ² / 230 g 144 cm ² / 67 g 100 cm ² / 50 g 46 cm ² / 40 g	S.M.P.S. (item 2.1 in the picture) Logic Board (item 2.2 in the picture) Inverter (item 2.3 in the picture) LCD Driver (item 2.4 in the picture) PCB Light Sensor (item 2.5 in the picture)
LCD	3.9 kg	18"

Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790


LCD screen FIMI MCL1802-IP10 / 9919-320-5072x



<div>Hazardous</div> <div></div> <div>To be Removed</div>	Substances:		Location
	Type	Quantity	
	Cd	0	
	Hg	21 mg max. (Mercury is present in the backlight lamps)	
	Pb	Lead is present in the soldering process of PCBs	
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	



Material		
Fe	2.3 kg	-
Al	0	-
Cu	0.1 kg	Cables
Plastics	1 kg	-
Boards (S ² > 10cm ²)	96 cm ² / 260 g	S.M.P.S. (item 2.1 in the picture)
	320 cm ² / 230 g	Logic Board (item 2.2 in the picture)
	144 cm ² / 67 g	Inverter (item 2.3 in the picture)
	100 cm ² / 50 g	LCD Driver (item 2.4 in the picture)
	46 cm ² / 40 g	PCB Light Sensor (item 2.5 in the picture)
LCD	3.9 kg	18"

LCD screen FIMI MML1801-GX / 9919-320-5122x


Hazardous  To be Removed	Substances:		Location
	Type	Quantity	
	Cd	0	
	Hg	21 mg max. (Mercury is present in the backlight lamps)	Next figure (1)
	Pb	Lead is present in the soldering process of PCBs	Next figure (2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	





Material		
Fe	6 kg	-
Al	0	-
Cu	0.11 kg	Cables
Plastics	1.5 kg	-
Boards (S² > 10cm²)	96 cm ² / 260 g 320 cm ² / 230 g 144 cm ² / 67 g 100 cm ² / 50 g 46 cm ² / 40 g	S.M.P.S. (item 2.1 in the figure) Logic Board (item 2.2 in the figure) Inverter (item 2.3 in the figure) LCD Driver (item 2.4 in the figure) PCB Light Sensor (item 2.5 in the figure)
LCD	3.97 kg	18"

Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790

LCD screen FIMI MML1802-GXR / 9919-320-5123x


Hazardous  To be Removed	Substances:		Location
	Type	Quantity	
	Cd	0	
	Hg	21 mg max. (Mercury is present in the backlight lamps)	Next figure (1)
	Pb	Lead is present in the soldering process of PCBs	Next figure (2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	




Material		
Fe	2.3 kg	-
Al	0	-
Cu	0.1 kg	Cables
Plastics	1 kg	-
Boards (S ² > 10cm ²)	96 cm ² / 260 g 320 cm ² / 230 g 144 cm ² / 67 g 100 cm ² / 50 g 46 cm ² / 40 g	S.M.P.S. (item 2.1 in the figure) Logic Board (item 2.2 in the figure) Inverter (item 2.3 in the figure) LCD Driver (item 2.4 in the figure) PCB Light Sensor (item 2.5 in the figure)
LCD	3.9 kg	18"


Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790


LCD screen PHILIPS FIMI MML1822-GXR / 9919-320-5136x

Hazardous  To be Removed	Substances:		Location
	Type	Quantity	
	Cd	0	
	Hg	21 mg max. (Mercury is present in the backlight lamps: 3.5 mg x 6 lamps)	
	Pb	0	
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	





LCD screen PHILIPS FIMI CML1812-GXR / 9919-320-5147x

Hazardous  To be Removed	Substances:		Location
	Type	Quantity	
	Cd	0	
	Hg	21 mg max. (Mercury is present in the backlight lamps: 3.5 mg x 6 lamps)	
	Pb	0	
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	




CRT screen FIMI TG17FM / 9896-010-0278x | PAGE 1 of 2

Hazardous  To be Removed	Substances		Location
	Type	Quantity	
	Cd	0	
	Hg	0	
	Pb	- Lead is present in the CRT glass - Lead is present in the soldering of PCBs	Next figure (item 1) Next figure (item 2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	
Special attention 	Item		Location
	<p>When handling or disposing of a CRT, you must take steps to avoid creating an implosion hazard for you or your trash removal service. The most simple and safe method to make the tube safe is to identify the small sealed glass nib at the far back of the tube (this may be obscured by the electrical connector) and then (while wearing safety glasses and gloves) filing a small nick across this and then to break it off using a pair of pliers. A loud sucking sound will be heard as the air enters the tube, releasing the vacuum. One must be very cautious not to break the neck of the tube when it is evacuated since there is no plastic coating preventing shattering of the glass. High vacuum and high voltage can be dangerous.</p>		Next figure (3)

CRT screen FIMI TG17FM / 9896-010-0278x | PAGE 2 of 2



Printed copies are uncontrolled




Material (kg)		
Fe	2.5 kg	-
Al	2.0 kg	-
Cu	1.3 kg	Cables
Plastics	4.4 kg	-
Boards ($S^2 > 10\text{cm}^2$)	cm ² 700 / 1860 g cm ² 650 / 1060 g cm ² 450 / 660 g cm ² 90 / 180 g cm ² 180 / 160 g cm ² 30 / 30 g	Power Supply (item 2.1 in the picture) Deflection Circuits (item 2.2 in the picture) Video Logic Board (item 2.3 in the picture) Magnetometer (item 2.4 in the picture) CRT Board (not visible in the picture) Keyboard (item 2.5 in the picture)
CRT	10.6 kg	21"

Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790

CRT screen FIMI FE17B / 9896-010-0296x | PAGE 1 of 2

Hazardous  To be Removed	Substances		Location
	Type	Quantity	
	Cd	0	
	Hg	0	
	Pb	- Lead is present in the CRT glass - Lead is present in the soldering of PCBs	Next figure (item 1) Next figure (item 2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	
Special attention 	Item		Location
	<p>When handling or disposing of a CRT, you must take steps to avoid creating an implosion hazard for you or your trash removal service. The most simple and safe method to make the tube safe is to identify the small sealed glass nib at the far back of the tube (this may be obscured by the electrical connector) and then (while wearing safety glasses and gloves) filing a small nick across this and then to break it off using a pair of pliers. A loud sucking sound will be heard as the air enters the tube, releasing the vacuum. One must be very cautious not to break the neck of the tube when it is evacuated since there is no plastic coating preventing shattering of the glass. High vacuum and high voltage can be dangerous.</p>		Next figure (3)

CRT screen FIMI FE17B / 9896-010-0296x | PAGE 2 of 2





Material (kg)		
Fe	6.9 kg	-
Al	1.0 kg	-
Cu	1.0 kg	Cables
Plastics	0.35 kg	-
Boards ($S^2 > 10\text{cm}^2$)	cm^2 77 / 80 g cm^2 550 / 360 g cm^2 100 / 280 g cm^2 788 / 1720 g cm^2 45 / 50 g	Raster Correction (item 2.1 in the picture) Video + CRT Board (item 2.2 in the picture) Mains Harmonic Reduction (item 2.3 in the picture) Mother Board (item 2.4 in the picture) Keyboard (item 2.5 in the picture)
CRT	7.7 kg	17"


Title: Recycling passport OD Eleva All Digital 708-028

DocID: XDR054-090790

CRT screen FIMI FE20B / 9896-010-0298x | PAGE 1 of 2



Hazardous  To be Removed	Substances		Location
	Type	Quantity	
	Cd	0	
	Hg	0	
	Pb	- Lead is present in the CRT glass - Lead is present in the soldering of PCBs	Next figure (item 1) Next figure (item 2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	
Special attention 	Item		Location
	When handling or disposing of a CRT, you must take steps to avoid creating an implosion hazard for you or your trash removal service. The most simple and safe method to make the tube safe is to identify the small sealed glass nib at the far back of the tube (this may be obscured by the electrical connector) and then (while wearing safety glasses and gloves) filing a small nick across this and then to break it off using a pair of pliers. A loud sucking sound will be heard as the air enters the tube, releasing the vacuum. One must be very cautious not to break the neck of the tube when it is evacuated since there is no plastic coating preventing shattering of the glass. High vacuum and high voltage can be dangerous.		Next figure (3)

CRT screen FIMI FE20B / 9896-010-0298x | PAGE 2 of 2




Material (kg)		
Fe	7.3 kg	-
Al	1.3 kg	-
Cu	1.1 kg	Cables
Plastics	0.5 kg	-
Boards ($S^2 > 10\text{cm}^2$)	cm ² 77 / 80 g cm ² 550 / 360 g cm ² 100 / 280 g cm ² 788 / 1720 g cm ² 45 / 50 g	Raster Correction (item 2.1 in the picture) Video + CRT Board (item 2.2 in the picture) Mains Harmonic Reduction (item 2.3 in the picture) Mother Board (item 2.4 in the picture) Keyboard (item 2.5 in the picture)
CRT	13.2 kg	20"

CRT screen FIMI TG21CM / 9896-010-0277x PAGE 1 of 2

Hazardous  To be Removed	Substances		Location
	Type	Quantity	
	Cd	0	
	Hg	0	
	Pb	- Lead is present in the CRT glass - Lead is present in the soldering of PCBs	Next figure (item 1) Next figure (item 2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	
Special attention 	Item		Location
	When handling or disposing of a CRT, you must take steps to avoid creating an implosion hazard for you or your trash removal service. The most simple and safe method to make the tube safe is to identify the small sealed glass nib at the far back of the tube (this may be obscured by the electrical connector) and then (while wearing safety glasses and gloves) filing a small nick across this and then to break it off using a pair of pliers. A loud sucking sound will be heard as the air enters the tube, releasing the vacuum. One must be very cautious not to break the neck of the tube when it is evacuated since there is no plastic coating preventing shattering of the glass. High vacuum and high voltage can be dangerous.		Next figure (3)

CRT screen FIMI TG21CM / 9896-010-0277x | PAGE 2 of 2



Printed copies are uncontrolled



Material (kg)		
Fe	4.5 kg	-
Al	0.2 kg	-
Cu	1.5 kg	Cables
Plastics	4.5 kg	-
Boards (S ² > 10cm ²)	cm ² 700 / 1920 g cm ² 650 / 1080 g cm ² 500 / 780 g cm ² 52 / 170 g cm ² 30 / 150 g cm ² 30 / 30 g	Power Supply (item 2.1 in the picture) Deflection Circuits (item 2.2 in the picture) Video Logic Board (item 2.3 in the picture) Magnetometer (item 2.4 in the picture) CRT Board (not visible in the picture) Keyboard (item 2.5 in the picture)
CRT	15 kg	21"


Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790

CRT screen FIMI MCD21CM / 9896-010-0287x | PAGE 1 of 2

Hazardous  To be Removed	Substances		Location
	Type	Quantity	
	Cd	0	
	Hg	0	
	Pb	- Lead is present in the CRT glass - Lead is present in the soldering of PCBs	Next figure (item 1) Next figure (item 2.x)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	
Special attention 	Item		Location
	<p>When handling or disposing of a CRT, you must take steps to avoid creating an implosion hazard for you or your trash removal service. The most simple and safe method to make the tube safe is to identify the small sealed glass nib at the far back of the tube (this may be obscured by the electrical connector) and then (while wearing safety glasses and gloves) filing a small nick across this and then to break it off using a pair of pliers. A loud sucking sound will be heard as the air enters the tube, releasing the vacuum. One must be very cautious not to break the neck of the tube when it is evacuated since there is no plastic coating preventing shattering of the glass. High vacuum and high voltage can be dangerous.</p>		Next figure (3)

CRT screen FIMI MCD21CM / 9896-010-0287x | PAGE 2 of 2

Printed copies are uncontrolled






Material (kg)		
Fe	4.0 kg	-
Al	2.5 kg	-
Cu	2.0 kg	Cables
Plastics	4.8 kg	-
Boards ($S^2 > 10\text{cm}^2$)	cm ² 700 / 1860 g cm ² 650 / 1060 g cm ² 450 / 660 g cm ² 90 / 180 g cm ² 180 / 160 g cm ² 30 / 30 g	Power Supply (item 2.1 in the picture) Deflection Circuits (item 2.2 in the picture) Video Logic Board (item 2.3 in the picture) Magnetometer (item 2.4 in the picture) CRT Board (not visible in the picture) Keyboard (item 2.5 in the picture)
CRT	16 kg	21"


Title: Recycling passport OD Eleva All Digital 708-028
 DocID: XDR054-090790

Cabinets:

FOLLOWING PAGES PROVIDE INFORMATION ON VARIOUS CABINETS POSSIBLY PRESENT IN THE SYSTEM.

DI2 I-CABINET FL-CPD-XTV(E), 0744-504-001 (9896 010 3308x and 9896 010 3309x)

Recycle Info	Items:	Location
Batteries  To be Removed	1x CR2032 3.0V Lithium coin cell of 2.8 gram	See next figure; slot BLA24 
Hazardous  To be Removed	Substances: BeCu (BerylliumCopper) contact springs Pb is present in the soldering process of some PCBs	See next figure








Battery on PCB in slot BLA24




BeCu contact springs



M-CABINET Eleva, 9896 001 4150x




Batteries  To be Removed	Type: 1x CR2032 3.0V Lithium coin cell	Location Next picture: 
	1x 3.5V Lithium battery	Next picture: 
Hazardous  To be Removed	Substances: BeCu (BerylliumCopper) contact springs	Location Next picture:  BeCu

Printed copies are uncontrolled




Title: Recycling passport OD Eleva All Digital 708-028
 DocID: XDR054-090790

Optimus TC 2T R/F 480V, 9890 000 6206x




Recycle Info	Items:	Location
Fluids / Gases 	<ul style="list-style-type: none"> Transformer oil, type: Shell Diala (This High Voltage-transformer oil contains no PCBs) 	See next figure (3)
Batteries 	Items: Lithium chrome cell, 3V Lithium chrome cell, 3V	Location See next figure (1) See next figure (2)
To be Removed		
Hazardous 	Substances: Pb is present in the soldering process of PCBs	See next figure
To be Removed		

steel, iron	iron, low alloy (<5%)	130,785	KG
	iron, high alloy (>5%)	2,334	KG
steel, iron		133,119	KG
nonferrous metals and alloys	aluminum, -alloy	12,274	KG
	copper, -alloy	5,005	KG
	zinc, -alloy	1,576	KG
nonferrous metals and alloys		18,855	KG
glass / ceramics	glass	0,288	KG
glass / ceramics		0,288	KG
plastics / organic substances	oil	48,06	KG
	thermoplastic	5,952	KG
	thermoset	4,245	KG
	elastomer	0,232	KG
plastics / organic substances		58,489	KG
standard parts	other electronic powered devices	45,92	KG
	printed circuit boards	11,584	KG
	cables	6,733	KG
	power supplies / transformers	6,184	KG
	mounting parts, attaching part	2,03	KG
	fans	1,804	KG
	motors, pumps	0,52	KG
standard parts		74,775	KG
		285,526	KG

Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790




Velara 2T GFD 480V, 9890 000 6209x| Velara 2T GCF 480V, 9890 000 7030x


Recycle Info	Items:	Location
Fluids / Gases 	<ul style="list-style-type: none"> Transformer oil, type: Shell Diala (This High Voltage-transformer oil contains no PCBs) 	See next figure (1) & (4)
Batteries 	Items: Lithium chrome cell, 3V Lithium chrome cell, 3V	Location See next figure (2) See next figure (3)
Hazardous 	Substances: Pb is present in the soldering process of PCBs	See next figure


steel, iron	iron, low alloy (<5%)	136,263	KG
	iron, high alloy (>5%)	5,674	KG
steel, iron		141,937	KG
nonferrous metals and alloys	aluminum, -alloy	12,655	KG
	copper, -alloy	5,075	KG
	zinc, -alloy	1,652	KG
nonferrous metals and alloys		19,382	KG
glass / ceramics	glass	0,288	KG
glass / ceramics		0,288	KG
plastics / organic substances	oil	64,06	KG
	thermoplastic	6,524	KG
	thermoset	4,505	KG
	elastomer	0,247	KG
plastics / organic substances		75,336	KG
standard parts	other electronic powered devices	45,92	KG
	printed circuit boards	11,835	KG
	cables	9,173	KG
	power supplies / transformers	6,184	KG
	mounting parts, attaching part	2,17	KG
	fans	1,804	KG
	motors, pumps	0,52	KG
	computers and accessories	0,218	KG
standard parts		77,824	KG
		314,767	KG

Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790



P-CABINET OD System, 9896-001-4172x

Hazardous  To be Removed	Substances:	Location
	BeCu (BerylliumCopper) contact springs	See next figure
	Pb is present in the soldering process of a PCB	See next figure



Title: Recycling passport OD Eleva All Digital 708-028
DocID: XDR054-090790